

# COVID-19 HOME TESTS

The following are authorized at-home over the counter (OTC) COVID-19 diagnostic tests under Emergency Use Authorization (EUA) by the FDA. All tests on this list are performed via nasal swab for self-testing.

Device	Company	Age	Time to Result	Clinical performance versus RT-PCR	Notes
<b>ANTIGEN TESTS</b>					
<b>CareStart COVID-19 Antigen Home Test</b>	Access Bio	≥2yrs	10mins	Positive agreement 87% Negative agreement 98%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>QuickVue At-Home OTC COVID-19 Test</b>	Quidel	≥2yrs	10mins	Positive agreement 83.5% Negative agreement 99.2%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 6 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>BinaxNOW COVID-19 Ag Card Home Test</b>	Abbott Diagnostics	≥2yrs	15mins	Positive agreement 84% Negative agreement 98.3%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> <li>Telehealth proctor supervision required.</li> <li>Smartphone or computer required.</li> </ul>
<b>BinaxNOW COVID-19 Antigen Self Test</b>	Abbott Diagnostics	≥2yrs	15mins	Positive agreement 84.6% Negative agreement 98.5%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>Flowflex COVID-19 Antigen Home Test</b>	ACON Laboratories	≥2yrs	15mins	Positive agreement 93% Negative agreement 100%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> <li>Alternative name: On/Go One COVID-19 Antigen Home Test.</li> </ul>
<b>BD Veritor At-Home COVID-19 Test</b>	Becton, Dickinson and Company	≥2yrs	15mins	Positive agreement 84.6% Negative agreement 99.8%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> <li>Smartphone required.</li> </ul>
<b>Celltrion DiaTrust COVID-19 Ag Home Test</b>	Celltrion USA	≥14yrs	15mins	Positive agreement 86.7% Negative agreement 99.8%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> <li>Mid-turbinate nasal swab.</li> </ul>
<b>Ellume COVID-19 Home Test</b>	Ellume Limited	≥2yrs	15mins	Positive agreement 95% Negative agreement 97%	<ul style="list-style-type: none"> <li>Use for individuals with or without symptoms.</li> <li>Mid-turbinate nasal swab.</li> <li>Smartphone required.</li> </ul>
<b>iHealth COVID-19 Antigen Rapid Test</b>	iHealth Labs	≥2yrs	15mins	Positive agreement 94.3% Negative agreement 98.1%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>MaximBio ClearDetect COVID-19 Antigen Home Test</b>	Maxim Biomedical	≥2yrs	15mins	Positive agreement 86.9% Negative agreement 98.9%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 5 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>CLINITEST Rapid COVID-19 Antigen Self-Test</b>	Siemens Healthineer	≥2yrs	15mins	Positive agreement 86.5% Negative agreement 99.3%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>SCoV-2 Ag Detect Rapid Self-Test</b>	InBios International	≥2yrs	20mins	Positive agreement 86.7% Negative agreement 100%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 5 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>COVID-19 At-Home Test</b>	SD Biosensor	≥2yrs	20mins	Positive agreement 95.3% Negative agreement 100%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 6 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>InteliSwab COVID-19 Rapid Test</b>	OraSure Technologies	≥2yrs	30mins	Positive agreement 85% Negative agreement 98%	<ul style="list-style-type: none"> <li>Use if symptoms began in the last 7 days, or if no symptoms, perform test 2 times over 3 days.</li> </ul>
<b>MOLECULAR TESTS</b>					
<b>Cue COVID-19 Test for Home and Over The Counter (OTC) Use</b>	Cue Health	≥2yrs	20mins	Positive agreement 97.4% Negative agreement 99.1%	<ul style="list-style-type: none"> <li>Use for individuals with or without symptoms.</li> <li>Smartphone and Cue Cartridge Reader required.</li> </ul>
<b>Lucira CHECK-IT COVID-19 Test Kit</b>	Lucira Health	≥2yrs	30mins	Positive agreement 92% Negative agreement 98%	<ul style="list-style-type: none"> <li>Use for individuals with or without symptoms.</li> </ul>
<b>Detect Covid-19 Test</b>	Detect	≥2yrs	60mins	Positive agreement 90.9% Negative agreement 97.5%	<ul style="list-style-type: none"> <li>Use if COVID-19 suspected, or if no symptoms, perform test 2 times over 3 days.</li> <li>Smartphone and Detect Hub required.</li> </ul>
<b>REFERENCES</b>					